

AIDA antibody

Catalog Number: 107850

Product name

AIDA antibody

Specificity

Human; other species not tested.

Antibody description

AIDA Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of AIDA recombinant protein (Accession Number: NM_022831). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107850(AIDA Antibody) at dilution of 1:600



IP Result of anti-AIDA (IP:Catalog No:107850, 3ug; Detection:Catalog No:107850 1:600) with HeLa cells lysate 4000ug.