

LARS2 antibody

Catalog Number: 112150

Product name

LARS2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LARS2 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in HepG2 cells, Jurkat cells, MCF7 cells. Observed molecular weight by Western-blot: 102 kDa

Preparation

This antibody was obtained by immunization of LARS2 recombinant protein (Accession Number: NM_015340). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112150(LARS2 antibody) at dilution of 1:400



IP Result of anti-LARS2 (IP:Catalog No:112150, 4ug; Detection:Catalog No:112150 1:500) with Jurkat cells lysate 2000ug.