

# Anti-OR10J1 antibody



Catalog Number: 176150

## Product name

Polyclonal

Anti-OR10J1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB, IF

Rabbit polyclonal antibody to OR10J1

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or10j1. aa range:241-290

IF: 1:100-1:500

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1:: Western Blot analysis of various cells using Olfactory receptor 10J1 Polyclonal Antibody