

# ARMC5 antibody



Catalog Number: 108150

## Product name

ARMC5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ARMC5 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, A2780 cells, HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of ARMC5 recombinant protein (Accession Number: BC014945). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108150(ARMC5 Antibody) at dilution of 1:300



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:108150(ARMC5 Antibody) at dilution of 1:300