

GFP tag antibody



Catalog Number: 110950

Product name

GFP tag antibody

Specificity

Tag

Antibody description

GFP tag Rabbit Polyclonal antibody. Positive IP detected in Transfected HEK-293 cells. Positive WB detected in recombinant protein. Positive IF detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of GFP tag recombinant protein (Accession Number: AVU73901). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP, IF

Dilutions


Recommended Dilution:

WB: 1:2000-1:10000

IP: 1:2000-1:20000

IF: 1:20-1:200


Validations



Western blot of eGFP-tagged fusion protein with anti-eGFP-tag (Catalog No:117319) at various

dilutions.

IP Result of anti-GFP tag (IP:Catalog No:117319, 4ug; Detection:Catalog No:117319 1:5000) with Transfected HEK-293 cells lysate 300ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed Transfected HEK-293 cells using Catalog

GFP tag antibody



Catalog Number: 110950

No:117319(GFP tag Antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)