

Human CSNK2A2 (His & GST Tag) recombinant protein



Catalog Number: 503950

General Information

Protein Construction

A DNA sequence encoding the human CSNK2A2 (NP_001887.1) (Met 1-Arg 350) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

The specific activity was determined to be 30 nmol/min/mg using casein as substrate.

Purity

> 92 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human CSNK2A2/GST chimera consists of 587 amino acids and has a calculated molecular mass of 69 kDa. It migrates as an approximately 60 kDa band in SDS-PAGE under reducing conditions.

Formulation

Supplied as sterile 50mM Tris, 500mM NaCl, 0.5mM PMSF, 0.5mM TCEP, 10% gly, pH 8.5

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CSNK2A2 / CK2A2 Protein (His & GST Tag) SDS-PAGE