

# PLCL2 antibody



Catalog Number: 113950

---

## Product name

PLCL2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PLCL2 Rabbit Polyclonal antibody. Positive IF detected in HepG2.

## Preparation

This antibody was obtained by immunization of PLCL2 recombinant protein (Accession Number: XM\_017006025). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IF

## Dilutions

Recommended Dilution:

IF: 1:10-1:100

## Validations



Immunofluorescent analysis of HepG2 cells, using PLCL2 antibody Catalog No:113950 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).