

# ZNF146 antibody



Catalog Number: 116950

## Product name

ZNF146 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ZNF146 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse lung tissue. Positive IHC detected in human heart. Observed molecular weight by Western-blot: 20-25kd, 33-35kd

## Preparation

This antibody was obtained by immunization of ZNF146 recombinant protein (Accession Number: BC005154). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:116950(ZNF146 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human heart using Catalog No:116950(ZNF146 antibody) at dilution of 1:100 (under 10x lens)