

Anti-Peroxiredoxin 2/PRP Rabbit antibody



Catalog Number: 167950

Product name

Rabbit IgG

Anti-Peroxiredoxin 2/PRP Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/2000-1/10000

Peroxiredoxin 2/PRP Rabbit monoclonal antibody

ICC/IF: 1/50

Preparation

FC: 1/20-1/100

Antigen: A synthetic peptide of human Peroxiredoxin 2/PRP

IP: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Validations


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of Peroxiredoxin 2/PRP in HeLa,A549,C6 using Peroxiredoxin 2/PRP Rabbit mAb(1:1000 diluted)