

Human CXCL6/GCP2/SCYB6 (His tag) recombinant protein



Catalog Number: 518950

General Information

Gene Name Synonym

C-X-C Motif Chemokine 6, Chemokine Alpha 3, CKA-3, Granulocyte Chemotactic Protein 2, GCP-2, Small-Inducible Cytokine B6, CXCL6, GCP2, SCYB6

Protein Construction

A DNA sequence encoding the Human CXCL6/GCP2/SCYB6 Gly38-Asn114 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 5% Trehalose, pH 7.4.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CXCL6/GCP2/SCYB6 (His tag)
recombinant protein