

Anti-RRH antibody

Catalog Number: 176250

Product name

Polyclonal

Anti-RRH antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB, IHC-P

Rabbit polyclonal antibody to RRH


Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rrrh. aa range:201-250

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Storage**

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Fig1.; Western Blot analysis of Jurkat cells using Peropsin Polyclonal Antibody

Clonality