

ZNF828 antibody

Catalog Number: 117250

Product name

ZNF828 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF828 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells.

Preparation

This antibody was obtained by immunization of ZNF828 recombinant protein (Accession Number: NM_032436). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:500-1:5000

Validations



IP Result of anti-ZNF828 (IP:Catalog No:117250, 4ug; Detection:Catalog No:117250 1:1000) with Jurkat cells lysate 4000ug.