

# Anti-beta Catenin (5D6) Mouse antibody



Catalog Number: 168250

## Product name

Anti-beta Catenin (5D6) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

beta Catenin (5D6) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human  $\beta$ -Catenin protein fragments expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of  $\beta$ -Catenin in HCT116, C6 and 3T3 cell lysates using  $\beta$ -Catenin mouse mAb(dilution 1:1000). Predicted band size: 92kDa. Observed band size: 92kDa.