

Human SDCBP/MDA9/SYCL (His tag) recombinant protein



Catalog Number: 519250

General Information

Gene Name Synonym

Syntenin-1, Melanoma differentiation-associated protein 9, Pro-TGF-alpha cytoplasmic domain-interacting protein 18, Scaffold protein Pbp1, Syndecan-binding protein 1, SDCBP, MDA9, SYCL,

Protein Construction

A DNA sequence encoding the Human SDCBP/MDA9/SYCL Ser2-Val298 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SDCBP/MDA9/SYCL (His tag) recombinant protein