

DDX31 Polyclonal Antibody



Catalog Number: 166650

Product name

DDX31 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to DDX31

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 652-851 of human DDX31 (NP_073616.6).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - DDX31 Polyclonal Antibody

Western blot analysis of extracts of A-549 cells, using DDX31 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 3min.