

Anti-FEN1 Rabbit antibody



Catalog Number: 167650

Product name

Anti-FEN1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

FEN1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human FEN1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/50-1/200

Validations



Western blot detection of FEN1 in HeLa, MCF-7, U251, HL-60 using FEN1 Rabbit mAb(1:1000 diluted)