

# C12orf72 antibody



Catalog Number: 108650

## Product name

C12orf72 antibody

## Specificity

Human; other species not tested.

## Antibody description

C12orf72 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of C12orf72 recombinant protein (Accession Number: NM\_001135864). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108650(C12orf72 Antibody) at dilution of 1:700