

CLDN4 Polyclonal Antibody



Catalog Number: 161350

Product name

CLDN4 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CLDN4

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human CLDN4 (NP_001296.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000


IHC 1:50 - 1:100

Validations



Western blot - CLDN4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CLDN4 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry - CLDN4 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using CLDN4 antibody at dilution of 1:150 (40x lens).