

# CD70 Polyclonal Antibody



Catalog Number: 163350

---

## Product name

CD70 Polyclonal Antibody

## Specificity

Mouse,Rat

## Antibody description

Polyclonal antibody to CD70

## Preparation

Antigen: A synthetic peptide of human CD70.

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - CD70 Polyclonal Antibody

Immunofluorescence analysis of rat spleen using CD70 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.