

# PXMP3 antibody



Catalog Number: 114350

## Product name

PXMP3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PXMP3 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 35-40 kDa

## Preparation

This antibody was obtained by immunization of PXMP3 recombinant protein (Accession Number: NM\_001079867). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114350(PXMP3 Antibody) at dilution of 1:600