

# SLC9A11 antibody



Catalog Number: 115350

## Product name

SLC9A11 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SLC9A11 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells. Observed molecular weight by Western-blot: 135-140kd

## Preparation

This antibody was obtained by immunization of SLC9A11 recombinant protein (Accession Number: NM\_178527). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115350(SLC9A11 antibody) at dilution of 1:500