

Insulin antibody

Catalog Number: 107350

Product name

Insulin antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Insulin Mouse Monoclonal antibody. Positive IHC detected in human pancreas tissue, mouse pancreas tissue.

Preparation

This antibody was obtained by immunization of Insulin recombinant protein (Accession Number: NM_001185097). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:500-1:5000

Validations



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:107350(Insulin Antibody) at dilution of 1:5000 (under 10x lens).



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using Catalog No:107350(Insulin Antibody) at dilution of 1:5000 (under 40x lens).