

Human Fetuin-B / FETUB (His Tag) recombinant protein



Catalog Number: 500750

General Information

Gene Name Synonym

16G2; Fetuin-like protein IRL685; Gugu

Protein Construction

A DNA sequence encoding the human FETUB (NP_055190.2) (Met 1-Pro 382) was expressed, with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met 19

Molecular Mass

The recombinant human FETUB consists of 375

amino acids and predicts a molecular mass of 41.7 kDa. rhFETUB migrates as an approximately 55 kDa band in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Fetuin-B / FETUB Protein (His Tag) SDS-PAGE