

RNF2 antibody

Catalog Number: 114750

Product name

RNF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNF2 Rabbit Polyclonal antibody. Positive IP detected in K-562 cells. Positive WB detected in human brain tissue, K-562 cells, NIH/3T3 cells, U-937 cells. Observed molecular weight by Western-blot: 38 kDa

Preparation

This antibody was obtained by immunization of RNF2 recombinant protein (Accession Number: XM_011509851). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114750(RNF2 antibody) at dilution of 1:300



IP Result of anti-RNF2 (IP:Catalog No:114750, 4ug; Detection:Catalog No:114750 1:800) with K-562 cells lysate 1200ug.