

# Anti-Phospho-eNOS (Ser1179) Rabbit antibody



Catalog Number: 166750

---

## Product name

Polyclonal

Anti-Phospho-eNOS (Ser1179) Rabbit antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat

## Applications

## Antibody description

IHC-P, WB

Phospho-eNOS (Ser1179) Rabbit polyclonal antibody

## Dilutions

## Preparation

IHC-p: 1/50-1/200

Antigen: Synthesized peptide derived from human eNOS (Phospho-Ser1179)

WB: 1/500-1/2000

## Formulation


PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality



Immunofluorescent analysis of Rat brain tissue, using eNOS Antibody .