

# ACOX2 antibody

Catalog Number: 107750

## Product name

ACOX2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ACOX2 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive IF detected in HepG2.

## Preparation

This antibody was obtained by immunization of ACOX2 recombinant protein (Accession Number: XM\_005265505). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF


## Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


## Validations



Immunohistochemical of paraffin-embedded human liver using Catalog No:107750(ACOX2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:107750(ACOX2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using ACOX2 antibody Catalog No:107750 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).