

Human CA5B (His tag) recombinant protein



Catalog Number: 518750

General Information

Gene Name Synonym

Carbonic Anhydrase 5B Mitochondrial, Carbonate Dehydratase VB, Carbonic Anhydrase VB, CA-VB, CA5B

Protein Construction

A DNA sequence encoding the Human CA5B Cys34-Pro317 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human CA5B (His tag) recombinant protein