

DCK antibody



Catalog Number: 109750

Product name

DCK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCK Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, Jurkat cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of DCK recombinant protein (Accession Number: NM_000788). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109750(DCK antibody) at dilution of 1:1000