

# ZIP4 antibody



Catalog Number: 117060

## Product name

ZIP4 antibody

## Clonality

Polyclonal

## Immunogen

[Human SLC39A4 Recombinant protein \(His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

Human; other species not tested.

## Applications

ELISA, WB

## Antibody description

ZIP4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, human jejunum tissue, MCF7 cells. Observed molecular weight by Western-blot: 60 kDa

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of ZIP4 recombinant protein (Accession Number: NM\_130849). Purification method: Antigen affinity purified.


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117060(SLC39A4 antibody) at dilution of 1:500