

ZIP4 antibody



Catalog Number: 117060

Product name

ZIP4 antibody

Immunogen

[Human SLC39A4 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZIP4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, COLO 320 cells, human jejunum tissue, MCF7 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of ZIP4 recombinant protein (Accession Number: NM_130849). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117060(SLC39A4 antibody) at dilution of 1:500