

# Anti-ACER2 antibody



Catalog Number: 176460

## Product name

Anti-ACER2 antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse

## Applications

WB, IHC-P

## Antibody description

Rabbit polyclonal antibody to ACER2

## Dilutions

WB: 1:500-1:2000

## Preparation


This antigen of this antibody was synthesized peptide derived from asah3l . at aa range: 50-130

IHC-P: 1:100-1:300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations



## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

Fig1:; Western blot analysis of COS7 using ASAH3L antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000