

Anti-ACER2 antibody



Catalog Number: 176460

Product name

Anti-ACER2 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ACER2

Preparation

This antigen of this antibody was synthesized peptide derived from asah3l . at aa range: 50-130

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western blot analysis of COS7 using ASAH3L antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000