

# Anti-Rab9 Rabbit antibody



Catalog Number: 167460

## Product name

Anti-Rab9 Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/50

## Specificity

Human, Mouse, Rat

ICC/IF: 1/50


## Antibody description

Rab9 Rabbit polyclonal antibody

## Validations

## Preparation

Antigen: A synthetic peptide of human Rab9



## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Western blot detection of Rab9 in 3T3,HeLa cell lysates using Rab9 Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF

## Dilutions



Immunohistochemistry of Rab9 in paraffin-embedded Human colon cancer tissue using Rab9 Rabbit pAb at dilution 1/50