

Human IL33/C9orf26/IL1F11/NFHEV recombinant protein



Catalog Number: 518460

General Information

Gene Name Synonym

Interleukin-33, IL-33, Interleukin-1 Family Member 11, IL-1F11, Nuclear Factor From High Endothelial Venules, NF-HEV, IL33, C9orf26, IL1F11, NFHEV

Protein Construction

A DNA sequence encoding the Human IL33/C9orf26/IL1F11/NFHEV Ser112-Thr270 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human IL33/C9orf26/IL1F11/NFHEV recombinant protein