

Human ECD / Ecdysoneless homolog (His & GST Tag) recombinant protein



Catalog Number: 500860

General Information

Gene Name Synonym

Protein ecdysoneless homolog; Suppressor of GCR2

Protein Construction

A DNA sequence encoding the human ECD (O95905-1) (Met1-Thr621) was expressed with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human ECD/GST chimera consists of 881 amino acids and has a calculated molecular mass of 100.6 kDa. The recombinant protein migrates as an approximately 88-108 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human ECD / Ecdysoneless homolog Protein (His & GST Tag) SDS-PAGE