

GAMT antibody



Catalog Number: 110860

Product name

GAMT antibody

Specificity

Human, Mouse, Rat, Zebrafish; other species not tested.

Antibody description

GAMT Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, human brain tissue, human heart tissue, mouse brain tissue, rat liver tissue. Positive IHC detected in human prostate cancer tissue, human liver tissue. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of GAMT recombinant protein (Accession Number: NM_000156). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110860(GAMT antibody) at dilution of 1:6000



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:110860(GAMT antibody) at dilution of 1:10 (under 25x lens)