

Integrin alpha-5 antibody



Catalog Number: 111860

Product name

Integrin alpha-5 antibody

Specificity

Human, Mouse, Monkey; other species not tested.

Antibody description

Integrin alpha-5 Rabbit Polyclonal antibody.
Positive WB detected in 3T3-L1 cells, A549 cells, COS-7 cells, HeLa cells, human placenta tissue.
Observed molecular weight by Western-blot:
120-135 kDa

Preparation

This antibody was obtained by immunization of Integrin alpha-5 recombinant protein (Accession Number: NM_002205). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:111860(ITGA5 antibody) at dilution of 1:1000