

TAZ antibody

Catalog Number: 115860

Product name

TAZ antibody

WB: 1:200-1:1000

IP: 1:200-1:1000

Specificity


Human; other species not tested.

IHC: 1:20-1:200

Antibody description

TAZ Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human heart tissue, human thyroid cancer tissue. Positive IP detected in MDA-MB-453s cells. Positive WB detected in MDA-MB-453s cells, LN229 cells. Observed molecular weight by Western-blot: 50 kDa

Validations




Preparation

This antibody was obtained by immunization of TAZ recombinant protein (Accession Number: NM_001168280). Purification method: Antigen affinity purified.

MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:115860(WWTR1,TAZ Antibody) at dilution of

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP


Dilutions

Recommended Dilution:

WB result of WWTR1 antibody (TAZ antibody, Catalog No:115860, 1:1000) with LN229 cell lysate by Dr. Xiaobing He. (Lane 1, LN229 (GBM cell line) cell lysate. 90 ug; Lane 2, LN229 growth inhibition treatment. 90 ug)



TAZ antibody

Catalog Number: 115860



Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:115860 (TAZ Antibody) at dilution of 1:50 (under 40x lens)

Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:115860 (TAZ Antibody) at dilution of 1:50 (under 10x lens)



IP Result of anti-WWTR1,TAZ (IP:Catalog No:115860, 4ug; Detection:Catalog No:115860 1:300) with MDA-MB-453s cells lysate 2800ug.