

Anti-MLH1 Rabbit antibody



Catalog Number: 168860

Product name

Anti-MLH1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

MLH1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MLH1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

1:1000-1:5000

Validations



Western blot detection of MLH1 in Rat Brain lysates using MLH1 Rabbit mAb(1:1000 diluted).Predicted band size:85kDa.Observed band size:85kDa.