

LBP antibody



Catalog Number: 112160

Product name

Polyclonal

LBP antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

LBP Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 53-60kd

Dilutions

Recommended Dilution:


Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of LBP recombinant protein (Accession Number: BC022256). Purification method: Antigen affinity purified.

Validations

Formulation



PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112160(LBP antibody) at dilution of 1:1000

Clonality