

ZNF320 antibody



Catalog Number: 117160

Product name

ZNF320 antibody

Immunogen

[Human ZNF320 Recombinant protein \(His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

ZNF320 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, L02 cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of ZNF320 recombinant protein (Accession Number: XM_047438304). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:117160(ZNF320 Antibody) at dilution of 1:300