

ARNTL antibody

Catalog Number: 108160

Product name

ARNTL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ARNTL Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, 3T3-L1 cells, HT-1080 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 69-75 kDa

Preparation

This antibody was obtained by immunization of ARNTL recombinant protein (Accession Number: NM_001351804). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108160(ARNTL Antibody) at dilution of 1:600



IP Result of anti-ARNTL (IP:Catalog No:108160, 4ug; Detection:Catalog No:108160 1:300) with HeLa cells lysate 1600ug.