

# CDK3 antibody



Catalog Number: 109160

---

## Product name

CDK3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CDK3 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, Jurkat cells, mouse brain tissue. Observed molecular weight by Western-blot: 40-45kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001258). Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109160(CDK3 antibody) at dilution of 1:800