

# Rat Il2/Il-2 (His tag) recombinant protein



Catalog Number: 519160

## General Information

### Gene Name Synonym

Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin; IL2

### Protein Construction

A DNA sequence encoding the Rat Il2/Il-2 Ala21-Gln155 is expressed with a 6His tag at the C-terminus.

### Organism

Rat

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 5mM  $\text{NaH}_2\text{PO}_4$ , 5mM Citric acid, 150mM NaCl, pH4.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Rat Il2/Il-2 (His tag) recombinant protein