

Rat Il2/Il-2 (His tag) recombinant protein



Catalog Number: 519160

General Information

Gene Name Synonym

Interleukin-2; IL-2; T-cell growth factor; TCGF; Aldesleukin; IL2

Protein Construction

A DNA sequence encoding the Rat Il2/Il-2 Ala21-Gln155 is expressed with a 6His tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 5mM NaH_2PO_4 , 5mM Citric acid, 150mM NaCl, pH4.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Rat Il2/Il-2 (His tag) recombinant protein