

# HSF4 antibody



Catalog Number: 111560

## Product name

HSF4 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HSF4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue, mouse brain tissue, mouse lung tissue, rat brain tissue. Observed molecular weight by Western-blot: 62-66 kDa

## Preparation

This antibody was obtained by immunization of HSF4 recombinant protein (Accession Number: NM\_001538). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111560(HSF4 antibody) at dilution of 1:1000