

UGT3A2 antibody



Catalog Number: 116560

Product name

UGT3A2 antibody

Specificity

human; other species not tested.

Antibody description

UGT3A2 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 50-55kd

Preparation

This antibody was obtained by immunization of UGT3A2 recombinant protein (Accession Number: NM_174914). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116560(UGT3A2 antibody) at dilution of 1:500