

S100A3 antibody



Catalog Number: 114960

Product name

S100A3 antibody

Specificity

human, mouse; other species not tested.

Antibody description

S100A3 Rabbit Polyclonal antibody. Positive WB detected in A2780 cells, HeLa cells, mouse skin tissue. Observed molecular weight by Western-blot: 7 kDa, 12 kDa

Preparation

This antibody was obtained by immunization of S100A3 recombinant protein (Accession Number: NM_002960). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A2780 cells were subjected to SDS PAGE followed by western blot with Catalog No:114960(S100A3 antibody) at dilution of 1:300