

# ZNF215 antibody



Catalog Number: 116960

## Product name

ZNF215 antibody

## Immunogen

[Human ZNF215 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF215 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 50-60 kDa

## Preparation

This antibody was obtained by immunization of ZNF215 recombinant protein (Accession Number: NM\_001354854). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116960(ZNF215 Antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:116960(ZNF215 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)