

# Anti-Geminin Rabbit antibody



Catalog Number: 168960

## Product name

Anti-Geminin Rabbit antibody

## Specificity

Human

## Antibody description

Geminin Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Geminin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Geminin in Jurkat cell lysates using Geminin Rabbit mAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.